# S · E · Y · S · S · O

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**RoHS**  
Compliant

## PRODUCT CONTENTS

- 1 brush head
- 1 brush handle
- 1 travel case with brush head disinfection and USB-C power input
- 1 USB-C cable

## PLUGGING IN AND CHARGING

- The sonic toothbrush and travel case can be charged using the USB-C cable.
- The first full charging of the toothbrush takes about 3 hours and the full charge will last for up to 12 weeks (if using the Sensitive mode), which means approximately 190 brushing cycles.
- Each next charging of the toothbrush should take about 2–3 hours.
- The first full charging of the travel case takes about 2.5 hours and the full charge will last for about 200 minutes of operation (40 disinfection cycles).
- Each next charging of the travel case should take about 1.5–2.5 hours.
- During charging, the toothbrush can be disconnected from the USB-C cable for brushing the teeth.
- After two minutes of brushing, you may reconnect the toothbrush and continue charging.
- There is no risk of overcharging the battery.
- To extend the battery life, it is recommended that the toothbrush should be charged when the LCD display shows a low battery level.

## USING THE SEYSSO GO SONIC TOOTHBRUSH

Install the brush head on the brush handle. Next, apply an appropriate amount of toothpaste to the bristles of the brush. Before turning on the toothbrush, place its head in your mouth. Switch on the toothbrush. The toothbrush turns off automatically after the selected mode's time elapses, and it is ready to be used again. While brushing, the toothbrush indicates an instruction to move to the next section of the teeth every 30 seconds.

The brush can operate in one of three speeds depending on

the selected mode: Clean (70,000), Sensitive (70,000), White (74,000).

## CHANGING THE BRUSH MODE

The toothbrush power button also controls the brushing mode selection. This allows you to cycle through the mode and brush speed in succession: Clean, Sensitive and White. Use the power button to turn on the toothbrush, then repeatedly press the power button to cycle through the brushing modes. If the selected mode is left on for more than 5 seconds, pressing the power button again turns off the toothbrush.

The SEYSSO Superior sonic toothbrush is turned on for the first time in the Clean mode.

- To change this mode, press the power button again.
- Keep pressing the power button repeatedly to cycle to the desired mode.
- The currently selected mode is indicated by the LCD display on the toothbrush.
- After 5 seconds, the selected mode is stored in the memory, after which the toothbrush will be turned off if you press the power button again.
- When the toothbrush is turned on again, it will run the last selected and memorised mode.

## WARNING! READ BEFORE USE!

KEEP THIS MANUAL FOR FURTHER REFERENCE.

- For home use only.
- The toothbrush features a waterproof handle. It is safe and designed to be used in a bathroom. The appliance has the IPX7 protection rating. It means that it is resistant to accidental immersion in water no deeper than 1 m for no more than 30 minutes.
- Remember: The toothbrush is not suitable for operating underwater!
- Do not attempt any unauthorised repairs.
- Straighten the cable before using it.
- In the event of a failure, do not disassemble the product or replace any parts, as there is a risk of electric shock. Contact the SEYSSO service centre.
- Replace the power cable if damaged. Contact the SEYSSO service centre.
- If dropped, replace the toothbrush head before next use. Do this even if there is no evident damage.

**UV LED light:**

The travel case features a UV-C radiation emitting LED. Direct exposure to UV-C light may have an adverse effect on your health. When closed tightly, the travel case protects against exposure to UV-C light. Do not stare directly into the UV-C light source, as it is visible through the cover. The LED light source is not a replaceable part. Do not attempt to remove it.

## MAINTENANCE AND CLEANING

- After brushing your teeth, rinse the toothbrush under running water.
- Remove the brush head from the brush handle and rinse it under running water.
- Do not use water hotter than 40°C or any acidic cleaning solutions.
- Wipe the brush handle dry before recharging it.
- Make sure the toothbrush and the travel case are disconnected from the power source before cleaning. Use only a dry cloth to wipe the outer surfaces of the charger clean!

**!** NOTE: DO NOT ATTEMPT TO WASH THE TRAVEL CASE UNDER RUNNING WATER OR SUBMERGE IT IN WATER.

- Keep your travel case in a clean and safe place.
- Do not connect the toothbrush or travel case to any power source unless completely dry.

Follow the basic precautions for handling and operating electrical equipment, especially with children nearby.

## BRUSH HEAD INSTALLATION AND REMOVAL

Before using the toothbrush for the first time, apply a thin layer of grease onto the brush handle metal drive end. This will make it easier to install and remove the toothbrush heads.

It is recommended to replace the brush head every 3 months.

## TRAVEL CASE

The travel case features UV-C light disinfection. The travel case can be charged using the USB-C cable.

Procedure:

Open the top of the travel case and place the brush head

and brush handle in the case. The disinfection cycle will begin when the travel case is closed.

To recharge the travel case, connect it to a power USB output or to the USB charger.

Charge for 1.5–2.5 hours.

\* The USB charger is not included.

The travel case has a standard USB port and is compatible with the SEYSSO Go sonic toothbrush handle.

#### Warning!

- Do not store the travel case where it can fall down or come into direct contact with water.
- Do not throw the travel case into water or submerge it in water. If the travel case gets wet, disconnect it from the power supply.
- Do not use it while taking a bath.
- Do not reach for the travel case which has fallen into water or other liquid. Unplug the travel case before cleaning and make sure that the case is dry before connecting it to a power source.
- Do not connect the travel case to a power source if the UV-C LED is still on with the case cover open. UV light can be harmful to the eyes and skin.
- Use a dry cloth to clean the travel case. Do not clean the case when the UV-C LED is still hot.

### DISINFECTING THE SEYSSO BRUSH HEAD WITH UV-C LIGHT IN THE TRAVEL CASE

The SEYSSO Go toothbrush travel case features brush head disinfection. The travel case is compatible only with the SEYSSO GO series brush heads. The UV-C light emitted by the LED eliminates bacteria, viruses and fungi present on the brush head fibres.

NOTE: Direct exposure to UV light is harmful to the eyes and skin.

- Only one toothbrush head can be disinfected at a time. The UV-C LED is located at the head tip.
- Place the brush head in the travel case. Make sure that it is inserted correctly, without its protective cap, and with the fibres towards the UV-C LED.

NOTE: To disinfect the brush head immediately after use, rinse it under running water and shake off excess water from the fibres.

- Close the travel case with the brush head inside, the disinfection cycle will start automatically and the UV-C LED

will come on.

**NOTE:** Do not use the travel case if it is impossible to close it fully as a result of damage.

- Full disinfection of the brush head takes 5 minutes. Opening the travel case earlier will end the disinfection cycle.

Disinfection for less than 5 minutes may be ineffective.

**CAUTION!** Do not use the travel case if the UV-C LED does not turn off when the case is opened. UV light can be harmful to the eyes and skin.

- Once the disinfection cycle is completed, the UV-C LED will turn off. If you want to repeat the disinfection process, open and close the case again.

**CAUTION!** Disinfection of brush heads with UV light may cause a gradual change in the fibre colour.

## SAFETY WARNING

To reduce the risk of electric shock:

- Always disconnect the product from the power source after use.
- Do not use it while taking a bath.
- Contact SEYSSO for detailed instructions of storage and use of the product and its included accessories.

To reduce the risk of fire, short-circuit or injury:

- Disconnect the charger from the power source when the charging process is completed.
- This product is not a toy. Take special caution if the toothbrush is used by or near a child.
- Use the product as intended. Use only the accessories recommended by the manufacturer.
- Do not use the charging cable if damaged. If the toothbrush does not work correctly after being dropped, left outdoors or submerged in water, return it to the original seller.
- Do not stretch the charger power cable, pinch it with doors or pass it over sharp edges. Do not turn on the toothbrush while it is charging. Keep the product components away from hot surfaces.
- Never use the product if you have problems concentrating or are feeling drowsy.
- Do not use the product outdoors or anywhere it can be exposed to aerosols or pure oxygen.
- Before you connect the product to a power outlet, verify its power specifications.

It is not recommended that children under the age of 3 years use the product. Children and people with limited sensory,

motor or mental abilities should use the product under the supervision of a guardian, or should be instructed in the correct and safe use of the product.

When brushing your teeth, do not apply excessive force or make excessively fast movements with the product. During the first days of using the toothbrush, the gums may bleed slightly. The bleeding should subside after a few days. If it persists or increases after this time, consult a dentist. For delicate gums and sensitive teeth, it is recommended to use the Sensitive mode.

## ENVIRONMENTAL PROTECTION

Do not dispose of the SEYSSO GO toothbrush together with household waste. Hand the product over to a waste electronic equipment disposal and recycling facility.

## ENVIRONMENTAL PROTECTION

PROBLEM	POSSIBLE CAUSE	SOLUTION
The toothbrush action seems to be weak.	The replaceable brush head might have been installed incorrectly on the handle.	Make sure that the brush head is properly installed.
		If the problem persists, replace the brush head with a new one.
After charging the toothbrush, the battery life is very short.	Incorrect connection of the toothbrush to the charger may result in ineffective charging.	Make sure that the toothbrush is properly connected to the charger.
		If the problem persists, the battery is probably reaching the end of its life. Contact the customer service department.
The charging LED does not come on when the toothbrush is connected to the charger.	If the toothbrush has not been used for over 6 months, it enters sleep mode.	After 6–20 minutes, the toothbrush should wake up from sleep mode and the battery charging should begin. This is normal.
		If the toothbrush still does not charge, contact the customer service department.
The toothbrush does not turn on.	Make sure if the power supply is sufficient.	Make sure that the toothbrush has been charged correctly.
		If the toothbrush does not turn on even after it has been recharged, contact the customer service department.
The UV-C LED does not come on/does not work.	The case cover is not completely closed.	Make sure that the case cover is completely closed.
	The UV-C LED is spent.	Contact the customer service department.

## ENVIRONMENTAL PROTECTION

PRODUCT DESIGNATION	SEYSSO GO TOOTHBRUSH
Voltage	3.7 V
Power	3 W
Switch type	Pushbutton
Mode display	LED indicators
Operating time	390 min maximum
Charging	Via USB-C port
Charging time	Approx. 2–3 hours
Battery	Li-ion, 800 mAh
Low battery indication	LED
Water ingress protection rating	IPX7

## SPECIFICATIONS OF THE SEYSSO GO TRAVEL CASE

PRODUCT DESIGNATION	SEYSSO GO
Voltage	3.7 V
Power	4 W
Light source type:	UV LED
Sterilising efficiency	99.9%
Disinfection cycle time	5 minutes
UV light wavelength	275 nm
Operating time	200 minutes
Charging	Via USB-C port
Charging time	Approx. 1.5–2.5 hours
Battery	Li-ion, 800 mAh
Low battery indication	LED

## TWO-YEAR WARRANTY

SEYSSO provides a warranty on this product for 2 years from the date of purchase. The warranty covers damage from manufacturing defects and incorrect workmanship of the product. Keep the proof of purchase for the duration of the warranty period. A warranty claim can be filed during the warranty period if the product has been used correctly and as intended.

### WARRANTY EXCLUSIONS:

Superficial damage resulting from normal use.

Damage caused by:

- a) accidental dropping or impact, misuse, abuse or commercial use, negligence or modification of the product,
- b) negligent maintenance or care,
- c) connection to a power source out of specifications,
- d) repairs or attempted repairs made by personnel not authorised by SEYSSO.

Any attempts to open or dismantle the product, its parts or accessories will render this warranty invalid, void and null. All warranty claims must be supported by evidence. Each warranty claim must be made within the guarantee period. The overlay of the brush handle is sensitive; it is recommended to handle it with care. Scratches or abrasions on the surface of the brush holder will be considered to be superficial damage from normal use. They are not covered by this warranty. The warranty applies to components that affect toothbrush performance.

SEYSSO HYPERTECH  
GO  
SONIC TOOTHBRUSH

WM Brands sp. z o.o.  
ul. Dziupli 12  
02-454 Warszawa, Poland  
Manufactured in P.R.C.

Model: SH.ST.020.BLK.1  
Voltage: 3.7 V  
Power: 3 W  
Lithium-ion battery

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**RoHS**  
Compliant

S·E·Y·S·S·O

REGISTER YOUR  
PRODUCT AT  
STREFASEYSSO.PL

Smile matters

WWW.SEYSSO.PL