S·E·Y·S·S·O

baby

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PLEASE READ THE INSTRUCTION LEAFLET CAREFULLY.

Children should be supervised while using the device.

The product is not toy. Store in place out of children's reach. Children shall not play, clean or repair the appliance. Toothbrush can be used by children and person with reduced physical, sensory or mental capabilities or lack experience or knowledge if they have been given supervision or instruction concerning the use of the appliance by person responsible for their safety.

Do not dip the device in water for too long. Do not alter the product any way other than battery brush head change. Do not place the product in high temperatures or expose to direct sunlight. Insert the battery with the correct polarity. Do not recharge the non-rechargeable batteries.

A rechargeable battery must be removed from appliance before being charged. Exhausted battery should be safely disposed of. Do not short-cut any circuits of the appliance.

HOW TO USE THE PENGUIN SONIC TOOTHBRUSH?

- Put a little amount of toothpaste on wet bristles.
- Press the button for light.
- Press the button one more time to start.
- Move the brush along the chosen arch for 30 seconds before moving on to the next one. You can make gentle sweeping movements. 30 seconds interval will let you know when to change the arch (of your baby's mouth). The Penguin Sonic Toothbrush will turn off after 2 minutes.
- Rinse the brush head with water and leave to dry to make it ready for next brushing time.

YOUR BABY'S SAFETY

This product is not a toy. Adult supervision is required at all time. The person responsible should perform the whole teeth cleaning. Keep out of children reach. Inspect carefully before each use. Remember to replace every brush head when it is damaged or worn. Do not chew. Do not submerge in water for a long time.

THE PRODUCT COMES WITH 2 BRUSH HEADS

The Newborn's Brush is for use with the baby between 0-18 months. The Tot's Brush is for use with the baby between 18-36 months. We recommended to use a new brush head after every 3 months for complete oral health.

CHANGING THE REPLACEABLE BRUSH HEAD

To detach: turn brush head anti-clockwise and pull upwise. To attach: align the new brush head and turn it clockwise. Note: notch on neck should align with inner notches of brush head to lock in place.

REPLACING THE BATTERY

- Turn the cap anti-counter to remove brush head.
 - Unscrew a screw with a screwdriver*.
 - Place the new battery AAA.
 - Tighten the screw.
 - Align line marker on battery cap with the line marker on brush handle.
 - Turn clockwise to replace the battery cap.
- *The screwdriver is not included in the set.

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