S · E · Y · S · S · O

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RoHS
Compliant

SET ELEMENTS

- 2 brush heads
- 1 handle with LCD screen
- 1 charging station
- 1 travel case with brush charging and brush head disinfection function, with USB port
- 1 USB cable

PLUGGING AND CHARGING

- Place the brush on the charger and insert the charger plug into a power supply socket.
- The LCD screen displays the battery charge level as a percentage.
- While charging, the value on the LCD screen percentage indicator increases over time.
- The first full charging process takes at least 8 hours and the battery holds the charge for up to 8 weeks (Sensitive mode), i.e. 110 brushing cycles.
- Each subsequent charging should take about 8 hours.
- While charging, the toothbrush may be removed from the charger for the time needed to brush your teeth.
- After brushing your teeth for 2 minutes, place the toothbrush on the charger again.
- There is no risk of overcharging the battery.
- To extend the battery life, it is recommended that the brush should be charged when the LCD screen indicates a low battery level.

SEYSSO SUPERIOR SONIC BRUSH SOFTWARE

Place the brush head on the brush holder. Next, apply an appropriate amount of toothpaste on the bristles of the brush. Before switching on the toothbrush, place its head in your mouth. Switch on the toothbrush. The appliance turns off automatically after the mode-controlled time elapses, and it is ready to be used again. While brushing your teeth, the interval occurring every 15 or 30 seconds informs you when to move to the next section of teeth being brushed.

The toothbrush operates at four speed settings, i.e. 70,000, 72,000 or 74,000 or 66,000 movements per minute, within the following modes: Daily (70 000), Sensitive (72 000), White (74

000), Polish (66 000) and Massage (70 000).

Operating time within individual brushing modes in minutes/seconds:

- Daily - 2 min (4 x 30 s);
- Sensitive- 2 min (4 x 30 s);
- White- 2,5 min (4 x 30 s + 2 x 15 s)
- Polish - 1 min (4 x 15 s);
- Massage- 1 min (4 x 15 s).

HOW TO CHANGE THE BRUSHING MODES?

The toothbrush activation button also controls the brushing mode selection. It facilitates switching the operation speeds sequentially between: Daily, Sensitive, White, Polish and Massage. Use the button to activate the appliance, then keep pressing it to move between brushing modes. If a mode is active for more than 5 seconds, the next use of the button switches off the appliance.

SEYSSO Superior sonic toothbrush is activated for the first time in the Daily. mode.

- In order to switch from Daily mode to Sensitive, mode, press the mode selection button again.
- To change an operation mode, press the function button.
- The current/selected operation mode is displayed on the LCD screen.
- After 5 seconds, the operating mode is stored in the memory, and the toothbrush is switched off after pressing the button again.
- When the toothbrush is activated again, it runs the last stored operation mode.

WARNINGS! READ BEFORE USING THE TOOTHBRUSH!

SAVE THIS MANUAL FOR FURTHER REFERENCE.

- For domestic use only.
- The toothbrush is equipped with a waterproof holder. It is safe and designed to be used in the bathroom ENVIRONMENT. The appliance comes with the IPX7 protection rating. It means that it is resistant to accidental submersion in water up to the depth of 1 metre for up to 30 minutes.
- Remember: the toothbrush is not suitable for operating underwater.
- Do not make any unauthorised repairs.

- Store the charger in a dry and safe place.
- Straighten the cable before using it.
- In the event of a failure, do not dismantle the appliance or replace any parts, as there is a risk of electric shock. Contact the SEYSSO service centre.
- Replace a damaged power cable. Contact the SEYSSO service centre.
- If dropped, replace the toothbrush head before its next use. This rule applies even if the damage is not visible.

Light UV-C:

The travel case comes with a light bulb. UV-C. Direct contact with UV-C light may have an adverse effect on the user's health. A tightly closed case protects against light UV-C, however, do not look at the light source, as it is visible through the cover. The light bulb is not a replaceable part, so do not remove it.

MAINTENANCE AND CLEANING

- After brushing your teeth, rinse the toothbrush under running water.
- Remove the brush head from the brush holder and rinse it under running water.
- Do not use water at temperatures exceeding 40°C or any acidic cleaning solutions.
- Wipe the brush holder, before placing it back on the charger.
- Before starting cleaning the charger, make sure it is unplugged from the power source. Use only a dry cloth to wipe the surface of the charger!

⚠ IMPORTANT: DO NOT WASH THE CHARGER UNDER RUNNING WATER, NOR SUBMERGE IT IN WATER!

- Store the charger at a clean and safe place.
- Do not plug the charger in until it is completely dry.

Take basic precautions with any electric appliance, especially when children are present.

BRUSH HEAD FITTING AND REMOVAL

Before using the toothbrush for the first time, apply a thin layer of lubricant onto the brush holder metal pin. This will make it easier to fit and remove the toothbrush heads.

We recommend changing brush these heads every 3 months.

LCD SCREEN

The brush comes with an LCD screen, which lights up when the brush is activated with the function button. The display shows:

- welcome screen;
- screen with a list of operation modes - to select a mode, press the function button "On/Off". Continue pressing the button to switch to another operation mode;
- screen with the currently selected operation mode, brushing time and battery level displayed while brushing teeth;
- brushing assessment screen displayed at the end of brushing;
- screen with the percentage level of battery charge displayed when the toothbrush is being charged.

TRAVEL CASE

The case also comes with a charger function. It facilitates charging the toothbrush via a USB cable.

Procedure:

- Open the top part of the case and place the brush holder inside.
- Plug the case to USB interface or to a USB charger plug*.
- Charge for 8 hours.

* the USB charger plug is not included in the set.

For the first 10 minutes of the toothbrush charging process, the brush head is disinfected. Disinfection is conducted regardless of the fact whether there is a brush head inside the case. After 10 minutes, the head disinfection ends, and the brush charging process continues.

The case comes with a standard USB port. It is compatible with the SEYSSO Superior sonic toothbrush holder.

Warnings!

- Do not store the case where it can fall down or come into direct contact with water.
- Do not drop or submerge the case in water. If the case gets wet, disconnect it from the power supply.
- Do not use it while taking a bath.
- Do not reach for a case which has fallen into water or other liquid. Unplug it before cleaning and make sure that the case is dry, before connecting it to a power source.
- Do not connect the travel case to a power supply if the UV light does not come on when the case is opened. UV light can be harmful to eyes and skin.
- Use a dry cloth to remove dirt from the case. Do not clean the case when the UV bulb is heated.

SEYSSO BRUSH HEAD DISINFECTION WITH UV LIGHT IN TRAVEL CASE

The SEYSSO Superior toothbrush travel case comes with a toothbrush head disinfection function. The case is compatible only with SEYSSO Superior series heads. The UV-A light emitted by the bulb eliminates bacteria, viruses and fungi present on the brush head fibres.

NOTE: Direct contact with UV light is harmful to eyes and skin.

- Only one toothbrush head can be disinfected at a time. The UV bulb is located at the head tip.
- Place the head in the case. Make sure that it is inserted correctly, without its protective cap, and with fibres directed towards the UV light.

NOTE: In order to disinfect the head just after use, rinse it under running water and shake excess water from its fibres.

- Close the case and connect it to a power supply, via a micro USB cable, to a laptop or power adapter.

NOTE! Do not use the case if it is impossible to close it fully as a result of damage.

- Full disinfection of the brush head takes 10 minutes. Opening the case ends the disinfection process. Disinfection lasting shorter than 10 minutes may be ineffective.

NOTE! Do not use the case if the UV light does not turn off when the case is opened. UV light can be harmful to eyes and skin.

- If there is a brush handle inside the case, it will charge simultaneously with the disinfection process. After the disinfection process is completed, the brush charging mode automatically remains active in the case.

NOTE! Disinfection of brush heads with UV light may cause a gradual change in the fibre colour.

WARNING - DANGER

In order to reduce the electric shock risk:

- Always disconnect the appliance from the power source after use.
- Do not use it while taking a bath.
- Contact SEYSSO for more detailed instructions related to storage and use of the appliance and the attached accessories.

In order to reduce the risk of catching fire, short-circuit or injury:

- Disconnect the charger from the power source when the charging process is completed.
- This product is not a toy. Take special caution when the toothbrush is being used by a child or in the presence of a child.
- Use it in accordance with its intended use. Use only accessories recommended by the manufacturer.
- Do not use the charger if the power cable or plug is faulty. If the toothbrush does not work correctly after being dropped, left in an open space, or submerged in water, send it back to the place of purchase.
- Do not stretch the charger power cable, jam it with a door or wrap it around sharp edges. Do not switch on the toothbrush while it is charging. Keep appliance components away from hot surfaces.
- Do not disconnect the charger from a power source by pulling on the cable. Use the plug for this purpose.
- Do not keep the plug with wet hands.
- Never use the product if you have problems concentrating and if you experience drowsiness.
- Do not use the appliance outdoors or in the places where it can be exposed to aerosols or pure oxygen.
- Plug the appliance into a mains socket. Read its specifications.

It is not recommended that children under the age of 3 years should use the device. Children and people with limited sensory, motor or mental abilities should use the appliance under the supervision of a guardian, or should be trained in its correct and safe use.

When brushing your teeth, do not use excessive force nor make excessively fast movements. During the first days of the toothbrush use, user's gums may bleed slightly. This bleeding should subside after a few days. If it persists or increases after this time, contact a dentist. For delicate gums and sensitive teeth, it is recommended to use the Sensitive mode.

ENVIRONMENTAL PROTECTION

Do not dispose of the SEYSSO toothbrush together with municipal waste. Hand it over to an electronic equipment disposal and recycling facility. The UV bulb contains mercury which can be harmful to the environment. Do not break it down or try to dismantle it.

TROUBLESHOOTING

PROBLEM	POSSIBLE CAUSE	SOLUTION
Toothbrush movement seems weaker.	The replacement brush head may have been incorrectly fitted on the brush holder.	Make sure that the brush head has been correctly fitted.
		If the problem reappears, replace the brush head with a new one.
After loading the toothbrush, the battery life is very short.	Incorrect positioning of the toothbrush on the charger can make the charging process ineffective.	Make sure that the toothbrush is placed vertically when charging.
		If the problem persists, battery wear is probably the cause. Contact the customer service department.
When the toothbrush is connected to the charging station, the charge level screen is not displayed.	If the toothbrush has not been used for over 6 months, it enters a sleep mode.	After 6 - 20 minutes, the toothbrush should awaken from the sleep mode and begin charging the battery. This is a standard behaviour.
		If the toothbrush still does not charge, contact the customer service department.
The toothbrush does not start.	Check if it is caused by insufficient supply power.	Check if the toothbrush has been charged correctly.
		If the toothbrush still does not start after charging, contact the customer service department.
The UV lamp does not light up/does not operate.	The case cover is not completely closed.	Make sure that the case cover is closed as required.
	The UV lamp is spent/worn.	Contact the customer service department.

SEYSSO SUPERIOR TOOTHBRUSH SPECIFICATIONS

PRODUCT NAME	SEYSSO SUPERIOR
Voltage	3.7 V
Power	3 W
Switch type	Button
Display	LCD
Information displayed	Operation modes, battery level, brushing assessment, brushing time
Operation time	Up to 230 minutes
Charging method	Inductive charging
Charging time	Approx. 5 - 8 h
Battery	Li-Ion 600 mAh
Battery-low alert	% indicator on LCD screen
Water resistance class	IPX7

TWO-YEAR GUARANTEE

SEYSSO provides a guarantee for the period of 2 years from the date of purchase. The guarantee covers damage related to manufacturing defects or incorrect workmanship of the product. Keep the proof of purchase for the duration of the guarantee period. A complaint can be filed during the guarantee period if the product has been used correctly and in line with its intended use.

THE GUARANTEE DOES NOT COVER:

Surface damage resulting from normal use.

Damage resulting from:

- a) accidental dropping or impact, misuse, abuse or commercial use, negligence or modification of the appliance,
- b) improper maintenance, lack of care of the product,
- c) connecting to an incorrect power source
- d) repairs or attempted repairs by a party not authorised by SEYSSO

Any attempts to open or dismantle the appliance, its parts or accessories will render the guarantee invalid. Guarantee claims must be supported by evidence. A claim must be made within the guarantee period. The material covering the brush holder is delicate, so it is recommended to exercise caution during use. Scratches or abrasions on the surface of the brush holder will be treated as surface damage resulting from normal use. Thus, they will not be covered by the guarantee. The guarantee applies to components that affect toothbrush operation.

SEYSSO HYPERTECH
SUPERIOR
SONIC TOOTHBRUSH

WM Brands sp. z o.o.
ul. Dziupli 12
02-454 Warszawa, Poland
Manufactured in P.R.C.

Model: SH.ST.010.GRP:1

Voltage: 3.7 V

Power: 3W

Lithium-ion battery

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www.seysso.pl

RoHS
Compliant

S·E·Y·S·S·O

REGISTER YOUR
PRODUCT AT
STREFASEYSSO.PL

Smile matters

WWW.SEYSSO.PL