S · E · Y · S · S · O

Instrukcja dotyczy modeli o numerach /
The instruction applies to models with the following numbers /
Návod sa týka modelov s číslami / Návod se týká modelů s čísly:
SCL.ST.005.STB.1, SCL.ST.005.MRO.1, SCL.ST.005.AGR.1,
SCL.ST.005.STB.2, SCL.ST.005.MRO.2, SCL.ST.005.AGR.2,
SCL.ST.005.PNK.1, SCL.ST.005.LAV.1, SCL.ST.005.BLU.1,
SC.ST.005.BLK.1, SC.ST.005.BLK.2, SC.ST.005.BLK.3GFTP,
SC.ST.005.BLK.3, SC.ST.005.WHT.1, SC.ST.005.WHT.2

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KIT COMPONENTS

- 1 toothbrush head
- 1 face washing tip in a case (models: SCL.ST.005.STB.2, SCL.ST.005.MRO.2, SCL.ST.005.AGR.2)
- 1 toothbrush case (models: SC.ST.005.BLK.3, SC.ST.005.WHT.2)
- 1 toothbrush handle
- Charging station with two slots for tips

CONNECTING AND CHARGING

- Place the toothbrush on the charger and connect the charger plug to the socket (socket specification: 100-240V- 50-60Hz).
- The charging indicator (Charge) light should turn on.
- A blinking (Charge) light indicates charging. If the light is steady, it means the toothbrush is fully charged.
- The first charge takes at least one day and lasts up to 6 weeks, or 84 brushing cycles.
- Each subsequent charge should take no more than one day.
- You can remove the toothbrush from the charger during charging to brush your teeth.
- After a two-minute brushing time, you can put the toothbrush back on the charger again.
- There is no risk of overcharging the battery.
- It is not recommended to discharge the battery completely; therefore, it is recommended to charge it every 4-5 weeks.

USING THE SONIC TOOTHBRUSH SEYSSO WITH A TOOTHBRUSHING TIP

Attach the tip to the handle. Then apply an appropriate amount of toothpaste to the bristles. Before turning on the toothbrush, place the tip in your mouth. Turn on the toothbrush. After two minutes, the device will turn off automatically and is ready for reuse. The toothbrush moves at two speeds, making 62,000 or 82,000 movements per minute in the modes: Clean (82,000), Soft (62,000), and Massage (62,000 – 82,000).

USING THE SONIC TOOTHBRUSH SEYSSO WITH A FACE WASHING AND MASSAGE TIP

Attach the tip to the handle. Then apply an appropriate amount of facial cleansing gel to the surface of the tip with soft, flexible silicone protrusions. Before turning on the toothbrush, hold it against your face so that it touches the skin surface. Turn on the toothbrush in Massage mode. Perform circular motions starting

from the center of the face towards the ears, jawline, and hairline. Avoid the eye and mouth areas. After two minutes, the device will automatically turn off and is ready for reuse. Rinse your face with lukewarm water. Then thoroughly wash the tip under running water. The other side of the tip is enriched with wide wavy protrusions with a convex structure designed for skin massage. You can apply cream or serum to your facial skin, and then restart the toothbrush in Massage mode. After two minutes, the device will automatically turn off and is ready for reuse. Then thoroughly wash the tip under running water.

HOW TO CHANGE WORK MODES

The toothbrush's power button also controls the selection of the working mode. With it, we can sequentially change the working speed between: Clean, Soft, and Massage. Use the button to turn on the device, then press it again to switch between the brush modes. If a mode is active for more than 5 seconds, the next press will turn off the device. The sonic toothbrush turns on in Clean mode.

- To switch from Clean mode to Soft mode, press the function selection button a second time.
- To switch to Massage mode, press it one more time.
- The mode indicator light shows which mode is currently being used.
- After 5 seconds, the working mode is memorized, and the next use of the button will turn off the toothbrush.

WARNINGS! READ BEFORE USE!

KEEP THESE INSTRUCTIONS For domestic use only.

- Place the toothbrush on the charger and connect the charger plug to the outlet (socket specification: 100-240V- 50-60Hz).
- The toothbrush is equipped with a waterproof handle. It is safe and designed for use in the bathroom. The toothbrush has an IPX7 protection rating. This means it is resistant to accidental immersion in water up to 1 meter deep for up to 30 minutes.
- Remember: the toothbrush is not designed to operate underwater!
- The product is not intended for self-repair.
- Store the charger in a dry and safe place.
- Straighten the cable before use.
- In case of failure, do not disassemble the device or replace any components, as there is a risk of electric shock. Contact

SEYSSO.

- A damaged power cord must be replaced with a new one. Contact SEYSSO.
- If dropped, replace the tip before next use. Even if the damage is not visible.

MAINTENANCE AND CLEANING

- After brushing your teeth or face, rinse the toothbrush under running water.
- Remove the tip from the handle and rinse it under running water.
- Wipe the handle before placing it back on the charger.
- Before cleaning the charger, make sure it is disconnected from the power supply. Use only a dry cloth to wipe the surface of the charger!

! **IMPORTANT: DO NOT WASH THE CHARGER UNDER RUNNING WATER OR IMMERSE IT IN WATER**

- Keep the charger in a clean and safe place.
- Do not connect the charger to power until it is completely dry.
- Keep the charger away from children.

APPLYING AND REMOVING THE TIP

Before the first use of the toothbrush, apply a thin layer of lubricant to the metal pin of the handle. This will facilitate the application and removal of tips.

We recommend replacing toothbrush tips every 3 months and face tips every 6 months.

WARNING ABOUT DANGER

To reduce the risk of electric shock:

- always disconnect the device from the power source after use.
- do not use while bathing.
- for more detailed instructions regarding storage and use of the device, as well as included accessories, refer to the manufacturer.
- To reduce the risk of ignition, short circuit, or injury:
- disconnect the charger from power after charging is complete.
- the product is not a toy. Exercise special caution if the toothbrush is used by or near children.
- use according to its purpose. Only use accessories recommended by the manufacturer.
- do not use the charger if the power cord or plug is damaged. If

the toothbrush does not work properly after being dropped, left outdoors, or immersed in water, return it to the place of purchase.

- do not stretch the charger's power cord, trap it in doors, wrap it around sharp edges. Do not turn on the toothbrush while charging. Keep device components away from heated surfaces.
- do not disconnect the charger from the power by pulling on the cable. Use the plug.
- do not handle the plug with wet hands.
- never use the product if you have problems concentrating or feel drowsy.
- do not use the device outdoors or in places where it can be exposed to aerosols or pure oxygen.
- connect the device to an outlet. Familiarize yourself with its specification.

It is not recommended to use the device by children under the age of 8. Children and persons with limited sensory, motor, or mental capacities should use the device under the supervision of a caregiver or should be properly trained on the correct and safe use of the device. Children should not play with the device. Unsupervised children should not perform cleaning and maintenance of the device.

TWO-YEAR WARRANTY

The manufacturer provides a warranty for a period of two years from the date of purchase. The warranty is valid in Poland and in the territory of official distributors. The warranty covers damages related to manufacturing defects or improper execution of the product. Keep the proof of purchase for the duration of the warranty. During the warranty period, you can file a complaint if the product was used properly and for its intended purpose.

THE WARRANTY DOES NOT COVER:

- 1 Surface damages resulting from normal use.
- 2 Damages resulting from:
 - a) accidental dropping or hitting, improper use, abuse, or commercial use, neglect, device modification,
 - b) improper maintenance, negligence towards the product,
 - c) connection to an incorrect power source,
 - d) repairs or repair attempts made by unauthorized persons.

Any attempts to open or disassemble the device, its parts, or accessories will void the warranty. Warranty claims must be

supported by proof. The claim must occur within the warranty period. The material covering the handle is delicate, so care should be taken during use. Scratches or wear on the surface of the handle will be treated as surface damages resulting from normal use and thus will not be covered by the warranty. The warranty applies to components that affect the functioning of the toothbrush. Proper care of your sonic toothbrush significantly increases its lifespan. Failure to follow recommendations may affect the loss of warranty.

This manual applies to models with numbers: SCL.ST.005.STB.1, SCL.ST.005.MRO.1, SCL.ST.005.AGR.1, SCL.ST.005.STB.2, SCL.ST.005.MRO.2, SCL.ST.005.AGR.2, SCL.ST.005.PNK.1, SCL.ST.005.LAV.1, SCL.ST.005.BLU.1, SC.ST.005.BLK.1, SC.ST.005.BLK.2, SC.ST.005.BLK.3GFTP, SC.ST.005.BLK.3, SC.ST.005.WHT.1, SC.ST.005.WHT.2

MANUFACTURER:
WM BRANDS SP. Z O.O.
UL. DZIUPIŁI 12 02-454
WARSAW, POLAND
Made in P.R.C.

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Charger input voltage: 100-240 V~, 50/60 Hz **Input current:** 30mA
Power: 1.7W **Water resistance (handle and toothbrush charger):** IPX7.
Protection against the effects of immersion in water up to 1 meter for 30 minutes. The toothbrushes have a built-in lithium-ion battery.

The product complies with the requirements of the so-called "New Approach" Directives of the European Union (EU).

The packaging is made of cardboard.
Suitable for recycling.

The manual is made of paper.
Suitable for recycling.

This device is marked in accordance with the Polish Act on Waste Electrical and Electronic Equipment and the WEEE Directive 2012/19/EU with the symbol of a crossed-out waste bin. This marking indicates that this equipment, after its service life, cannot be disposed of with other household waste. The user is obliged to hand it over to those conducting the collection of waste electrical and electronic equipment. Those conducting the collection, including local collection points, stores, and municipal units, create a proper system that allows the disposal of this equipment. Proper disposal of waste electrical and electronic equipment helps avoid harmful consequences to human health and the natural environment, resulting from the presence of hazardous substances and improper storage and processing of such equipment.

Symbol indicating Class II protection for electrical devices.

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